


DETERGENT - PRE-DISINFECTANT
MEDICAL DEVICE
CLASS IIb



UMONIUM³⁸ STERILY



✓ FORMULA FREE FROM

- CMR COMPONENT**
- ENDOCRINE DISRUPTOR

“ The number 1 in
sterilisation units ”



No need for enzyme or alkaline products to ensure sparkling instruments. UMONIUM³⁸ STERILY in the first step removes dirt, does not leave any yellow-brown deposits on the shaft of clamps and deep cleans instruments by soaking in an ultrasonic bath and/or by manual brushing. Safe for the operator, it quickly kills the population of micro-organisms. The solution is non-hazardous if it leaks, non-corrosive, biodegradable and does not contain any toxic vapours or substances likely to bind proteins. Equipment sparkles, waste pipes are clear of sticky waste and the wastewater is non-polluting.

** Carcinogenic, Mutagenic and Reprotoxic

 **laboratoire
HUCKERT'S**
INTERNATIONAL
mastered disinfection





✓ FORMULA FREE FROM

- CMR COMPONENT**
- ENDOCRINE DISRUPTOR

INDICATIONS

Cleaning and pre-disinfection for invasive medical devices before sterilisation. Prepare for disinfection.

COMPOSITION

N-benzyl-N-dodecyl-N, N-dimethyl-ammonium chloride/N-benzyl-N, N-dimethyl-N-tetradecyl-ammonium chloride. 22g/L.

REGULATORY STANDARDS

- Complies with ED 93/42/EEC relating to medical devices (annexe II, excl. §4)
- Quality and traceability guaranteed by ISO 9001:2015 and EN ISO 13485:2016


INSTRUCTIONS FOR USE

Use at 0.5% (25 mL/5 L of water*). Soak then scrub/use the ultrasounds to eliminate all dirt. Rinse thoroughly with clean water before disinfection.

PROPERTIES

- Medical device Class IIb
- Pre-disinfectant 0.5% cleaner
- Non-toxic, no dangerous vapours
- LD₅₀ (rats): 5840 mg/kg
- Non-corrosive: pH 7
- Not irritating to the skin (no transdermal toxicity)
- Aldehyde, peroxide and biguanide free, no other carcinogenic, mutagenic, reprotoxic components
- Non-coagulant
- Non-flammable
- 100% compatible on fibre optic, rubber, polycarbonate, acrylic compound, glass, Pyrex, PU, PVC, HDPE, PET, neoprene, latex, silicone, paints and stainless steel 410
- Effective, fast and efficient (dissolves blood)
- Effective in the presence of organic materials
- Faster acting than enzyme solutions
- No residue after rinsing
- Stability: 36 months date of manufacture
- Biodegradable

PROCEDURE FOR USE



Use at 0.5% (25 mL/5 L of water*)
→ Use of a bath by batch



Soak then scrub/use the ultrasounds to eliminate all dirt



Rinse thoroughly with clean water before disinfection



ITEM No.	NAME	BOX
PF 10320	1 litre dosing bottle	12 x 1 litre
PF 10322	5 L drum with pump	4 x 5 litres
PF 10324	125 L dosing bottle	24 x 125 mL
EMB 0013	ECONOMATIC pump complete kit - Dosage 0.5%	1 item

* Clean drinking water at room temperature

** Carcinogenic, Mutagenic and Reprotoxic

www.huckerts.net

FP00152 EN - Ed. 3.3 - 02/2019



BELGIUM
20 Avenue Lavoisier
1300 Wavre
Tel.: +32 (0)67 89 41 00
Fax: +32 (0)67 84 37 67

FRANCE
38 Avenue du Centre
78230 Le Pecq-sur-Seine
Tel.: +33 (0)1 3976 1505
Fax: +33 (0)1 3976 1115

G.D. LUXEMBOURG
19 Rue de l'Industrie
8069 Bertrange
Tel.: +352 26 39 42 60
Fax: +352 26 39 42 70

ITALY
20 Via Giulio Natta
05100 Terni
Tel.: +39 0744 80 03 87
Fax: +39 0744 81 73 39

SINGAPORE
51 Changi Business Park Central 2
The Signature, #04-05
Singapore 486066
Tel: +65 983 45 333



UMONIUM³⁸ is a registered trademark of Laboratoire Huckert's International