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

MADE IN BELGIUM

BACTERICIDAL, MYCOBACTERICIDAL, FUNGICIDAL, VIRUCIDAL

EN1040:2006, EN1275:2005, EN13727:2003, EN13624:2004, EN14348:2005, EN14476:2013, EN14561:2006, EN14562:2006, EN14563:2009

UMONIUM³⁸ NEUTRALIS

EN INDICATIONS: Cleaning and disinfection for invasive medical devices before sterilisation. Cleaning and disinfection for non-invasive medical devices and their surfaces. **INSTRUCTIONS: Cleaning/Pre-disinfection:** Use at 0.5% (25mL/5L of water*). Soak then scrub/use the ultrasounds or wet wipe to eliminate all dirt. Rinse thoroughly with water* or using a cloth dampened with water* if you are cleaning surfaces. **Disinfection:** Use at 2.5% (25 mL/L of water*) on clean devices. Apply by wet wiping or soaking, making sure to respect a minimum contact time of 10 minutes. If you are cleaning invasive medical devices: Rinse with water*. If you are cleaning non-invasive medical device surfaces: Leave to dry or rinse with water* (or using a cloth dampened with water*) if necessary. *Clean drinking water at room temperature. **PRECAUTIONS:** Pure product may cause irritation to eyes and skin. Avoid skin contact. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Do not discharge pure product into the environment. Consult specific directions and material safety data sheet.

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ET00V83 - Ed 4.0 - 09/2018

COMPOSITION UMONIUM³⁸ NEUTRALIS (SL):
N-benzyl-N-dodecyl-N,N-dimethyl-ammonium chloride / N-benzyl-N,N-dimethyl-N-tetradecyl-ammonium chloride. 99g/L.

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