

Tablet Hardness Tester

Precise measuring of up to 5 parameters Fast and intuitive touch screen operation Proven Dr. Schleuniger® hardness measuring Professional reporting including statistics



Swiss Quality High quality components paired with Swiss precision engineering provide for a long product life and minimize costly machine downtime. In order to ensure maximum reliability of your MultiTest 50, all instruments are tested according to strict in nal quality standards before they leave o **Hardness Testing** tablet breaking force) and EP (2.9.8, resistance to Quality matters.

MultiTest 50

The MultiTest 50 is a manual tablet testing system that offers highest precision, flexibility and exceptional user-friendliness. From basic hardness tests to complete measuring of up to five parameters, including professional reporting, physical testing has never been easier. Besides fast and efficient operation, the MultiTest 50 provides accurate and highly reliable test data, resulting in increased productivity, efficiency and ease of regulatory compliance. Its outstanding quality with proven Dr. Schleuniger® hardness measuring technology makes the MultiTest 50 the ideal choice for research laboratories, quality control departments and production facilities of pharmaceutical and nutraceutical companies.




Flexibility and Precision

Samples are placed inside the MultiTest 50 test area individually and by hand. Practically any imaginable design, including oblong, diamond, hexagon and other special shapes can be reliably positioned to measure hardness, diameter/length, width, thickness and weight. Special jaws and quick-change tongues allow testing of even the most unconventional tablet shapes.

The open and easily accessible test area allows for fast
The solid construction of the MultiTest 50 guarantees and reliable tablet orientation by the operator. To ensure exchanged in a matter of seconds. Even non-symmetrical shapes can be tested using a custom tongue.

reliable test results at all times. Dr. Schleuniger® hardcorrect alignment of special shapes without limiting the ness measuring technology complies with all current flexibility of the tester, the MultiTest 50 features quick- USP and EP requirements and offers consistently accuchange tongues with or without groove that can be easily rate results. Tablet hardness is measured using a precision S-beam load cell with a linearity of 99.95 % over the entire measuring range.



→ Testing of virtually all tablet shapes and sizes with the MultiTest 50









Fast, Faster, FastTest™

The MultiTest 50 is a manual multi-purpose instrument that optionally accumulates measurements from an external balance and thickness gauge in a single report. In FastTest™ configuration, operators can test tablets as quickly as with any automatic instrument.

of tablets, thickness is measured using the integrated operator to actuate testing for each tablet. precision gauge. With one hand turning the carousel

Testing larger sample sizes manually can be time-con- and the other controlling the thickness gauge, accurate suming and requires strict sequential testing by the measurements for the entire series are taken in very user to assure that all 5 measurements are correctly short time without compromising the order of test samassigned to individual tablets. In FastTest™ configurati- ples. Width, diameter/length and hardness testing is on, the operator transfers tablets after weighing to a executed afterwards by placing single tablets inside the carousel mounted on top of the MultiTest 50. Once test area. In automatic mode, individual tests are starweight measurements have been acquired for a series ted at programmable intervals without requiring the



Hardness testing \Rightarrow Manual tablet testing system \Rightarrow MultiTest 50

User-friendliness

Operating a manual tablet tester has never been easier. The intuitive icon-based touch screen user interface sets new standards for fast and user-friendly control. Three test modes easily retrieved for fast test setup. Up allow using the MultiTest 50 according to individual quality to 100 products including nominal requirements.

From executing simple tests to programming products with T1/T2/PL limits and methods - everything is available at the touch of a button. Intuitive icon-based navigation and logical menu structures with intelligent operator assistance provide for fast test setup, short training times and error-free data entries.

Different views that can be selected while test runs are executed allow the operator to display currently measured values, complete series of measurements and even graphical hardness curves. Professional reports are automatically printed at completion of a test run.

For paperless reporting and central data management in full compliance for each parameter. No programming with 21 CFR, Part 11 requirements, the MultiTest 50 can be networked with the TabStat™ DataCenter.





→ Fast and intuitive icon-based touch screen operation

"FullTest" mode

All product and method specifications can be pre-programmed and values and T1/T2/PL limits for all test parameters are conveniently stored directly on the MultiTest 50. On completion of each test run, a comprehensive report is printed automatically. The report shows all product and test parameters, calibration data and individually measured values with a statistical analysis.

"EasyTest" mode

For tests without T1/T2/PL limits, the EasyTest mode allows direct setting of nominal value and sample size of product or method specifications required. All reporting functions known from the "FullTest" mode remain available.

"SingleTest" mode


To quickly measure single tablets without printing a report, simply activate the parameter(s) you would like to test and press start.

Networking


All MultiTest 50 hardness testers can be networked with TabStat™ Data-Center software for 21 CFR, Part 11 compliant electronic data management. Product and method specifications are administrated centrally and simply loaded to individual testers via network when setting up a new test. Measured results and test reports are automatically stored in a central database. All reports are available anytime for printing and can be exported as PDF or in Microsoft Excel® compatible format.











ightarrow Professional data management and paperless reporting with TabStat[™] DataCenter

Technical Specifications

MultiTest 50, Tablet hardness tester

Weight	Measuring range	0.1 mg –100 g		
	Resolution	0.1 mg		
	Units of measure	mg/g		
	Type of external balance	Mettler Toledo, Sartorius		
Thickness	Measuring range	0 – 13 mm (0 – 0.51 inch)		
	Resolution	0.01 mm (0.00039 inch)		
	Units of measure	mm/inch		
Diameter / Length / Width	Measuring range	0-35 mm (0-1.37 inch); optional 65 mm (2.55 inch)		
	Resolution	0.01 mm (0.00039 inch)		
	Units of measure	mm/inch		
Hardness	Measuring range	Standard: up to 400 N Optional: 50 N / 500 N / 800 N		
	Accuracy	Better than +/- 1 N		
	Units of measure	N, Sc, Kp, programmable user-defined factor (based on N)		
Hardness measuring principle	Constant speed	0.05 – 5 mm/sec, programmable		
	Constant force (load rate)	10-200 N/sec, programmable		
	Constant power increase	Standard "6D" mode		
Interfaces		Ethernet, USB, RS232		
Power supply		100 – 230 V, 50 – 60 Hz		
CE Conformity		All MultiTest 50 versions fully comply with all CE and EMC equipment guidelines relative to mechanical and electrical safety and electromagnetic compatibility.		

		н	DH	TDH	WTDH	FastTest™
Configuration	Hardness	•	•	•	•	•
	Diameter / Length	_	•	•	•	•
	Width	_	•	•	•	•
	Thickness (internal)	_	_	•	•	•
	Thickness (external)	_	_	•	•	•
	Weight (external)	_	_	_	•	•
	FastTest™ configuration	_	_	0	0	•
Weight	(without packaging)	approx. 10 kg (22 lbs)	approx. 11 kg (24 lbs)			
Dimensions	Width	420 mm (16.5 inch)				
	Height	185 mm (7.3 inch)	185 mm (7.3 inch)	185 mm (7.3 inch)	185 mm (7.3 inch)	300 mm (11.8 inch)
	Depth	180 mm (7.1 inch)				
Software	TabStat™ DataCenter	0	0	0	0	0
Options	Function checking devices	0	0	0	0	0
	Calibration tools	0	0	0	0	0
	Custom tongues	0	0	0	0	0
	Special jaws	0	0	0	0	0

Technical specifications are subject to change without prior notice. Products illustrated in this brochure may include options or modifications not fitted as standard. No liability for errors and omissions.



Pharmatron AG | Uttigenstrasse 28 | CH-3600 Thun | Switzerland P +41 [0]33 227 50 00 | F +41 [0]33 227 50 01 | info@pharmatron.ch

Your local sales & service partner