## supplierfilterair.com



Each box of 25 vials is supplied with a Hanna certificate of quality. The reagents are traceable to NIST SRM® 930.

- Three measurement ranges to satisfy every need
  - As COD levels vary depending on the application and process measuring points, Hanna offers reagents to cover three separate ranges. Simply choose the best range for the application: low range: 0 to 150 mg/L O<sub>2</sub> medium range: 0 to 1500 mg/L O<sub>2</sub> high range: 0 to 15000 mg/L O<sub>2</sub>

## Accurate and repeatable measurements

- Hanna COD reagents have been developed in accordance with Standard Methods 5220D, USEPA 410.4 and ISO 15705:2002 methods.
- Pre-dosed vials
  - · Hanna vials contain approximately 3 mL of pre-dosed reagent. The operator just needs to add a small quantity of the sample.
- Ouick and accurate measurements
  - With pre-dosed vials, test preparation time is dramatically reduced. There is no time-consuming reagent preparation procedure or glassware cleaning.

## Safe reagents

· Hanna COD reagents are safe for operators and the environment. Vials and caps have been designed to avoid accidental reagent spills. Due to the pre-dosed reagents, the amount of chemicals and handling time is minimized.

COD Test	Range	Method	Reagent Code	
COD LR	0 to 150 mg/L	dichromate EPA‡ dichromate mercury-free°° dichromate ISO°	H193754A-25 H193754D-25 H193754F-25	
COD MR	0 to 1500 mg/L	dichromate EPA‡ dichromate mercury-free°° dichromate ISO°	HI93754B-25 HI93754E-25 HI93754G-25	
COD HR	0 to 15000 mg/L	dichromate	HI93754C-25	

## **COD** Standards

HI93754A

HI93754-11 500 ppm COD standard, 500 mL bottle HI93754-12 14000 ppm COD standard, 500 mL bottle

Notes: Method with chromium-sulfuric acid is officially recognized by EPA for wastewater analysis. This method is recommended for general purpose analysis with no chloride interference. Wethod follows the official method ISO 15705. COM RHS Do method is 0-1000 mg/L Meter can read higher.



www.hannainst.com



certified reagents

11

11.15