

FTC 120

Cooled Incubator

Designed for BOD incubation



supplierfilterair.com

ligh efficency with A+ energy rate

- Usable volume 109 L
- Uniform temperature by the use of forced ventilation (± 0,5°C)
- Constant temperature 20 °C
- Cooling system in A+ class energy rate
- Internal sockets to power devices
- 2-digit display of internal temperature
- Accepts 3 BOD Sensor System 6, or 2 BOD Sensor System 10
- Accepts 3 RESPIROMETRIC Sensor System 6



INSTRUMENT - CODE

FTC 120

230V/50-60Hz

F10300305



TECHNICAL DATA

DIGITAL DISPLAY	2-digit, 1 °C resolution
USABLE INTERNAL VOLUME	109 liters
ELECTRONIC THERMOREGULATION SYSTEM	Traditional
TEMPERATURE REGULATION	20 °C constant
INTERNAL TEMPERATURE STABILITY	± 0.5°C
INTERNAL TEMPERATURE HOMOGENEITY	± 0.5°C
NUMBER OF SHELVES INCLUDED	2
NUMBER OF INTERNAL SOCKETS	2
POWER	120 W
POWER SUPPLY	230 V / 50-60 Hz
WEIGHT	36 kg (79.4 lb)
DIMENSIONS (WXHXD)	540x912x550 mm (21.3x35.9x21.7 in)







We reserve the right to make technical alterations We do not assume liability for errors in printing, typing or transmission

VELP Official Partner



ITALY - HQ Via Stazione 16 20865 Usmate (MB) Italy Tel. +39 039 628811 velpitalia@velp.com

INDIA velpindia@velp.com

DESIGNED AND MANUFACTURED IN ITALY

155 Keyland Court, Bohemia NY 11716 - U.S. Tel. +1 631 573 6002 velpusa@velp.com

CHINA

Xinlong Rd Building 28, Lane 1333 Shanghai city - China Tel. +86 18017557329 velpchina@velp.com