

















Patient on the mattress can be moved parallelly to another bed or flat table

Specification:

- Construction
 - Made of steel pipe, square and steel sheet
- Hi-Low
 - Adjustable manually by crank
- Deck mattress
 - Made of multiplex
- Mattres
 - Made of acrylic laminated, can be moved
- Side guard

Made of aluminum alloy protective rails can be folded and hidden under the surface on no used

- Castor
 - Four castor Ø6", twin wheel, with central locking system
- Infuse stand
 - Available, Made of stainless steel pipe adjustable height with screw lock
- Painting
 - Epoxy powder coating anti bacterial grade
- Dimension (approx)
 - 1950 (L) x 630 (W) x 600 (H min) 820 (H max) mm

