

## Bed Series

**büma**<sup>®</sup>  
Medical Equipment



FK.01.04/1/007/2017



**ICU/ICCU Bed / 21021-EL**

**(Tempat Tidur / Ranjang Pasien ICU/ICCU)**

**KEMENKES RI AKD 20902410287**

e catalogue  
e-katalog.lhpp.go.id

Side Guard  
Can be Put Down



Low Position



High Position



Trendelenburg Position



Reverse  
Trendelenburg Position



### Specification :

#### ✓ Construction

Made of combination of plastic and steel sheet, square pipes steel and pipes steel. Assembly with a system of MIG, TIG and SPOT welding process.

#### ✓ Part

##### 1. Remote control :

Consist of 3 parts (sideguard on the right and left side, and on the outer foot panel.)

##### 2. Mattress platform

Made of plastic and steel, painting with powder coating and emboss ventilation holes, easy to clean, completed with holes for infuse stand.

##### 3. Head & Foot Panel

Polypropylene head & foot board, detachable & easy to clean.

##### 4. Castor

Four castor, swivel 360° with central locking system.

##### 5. Side guard

Made of plastic, 4 (four) separated section, automatic lock with gas spring.

#### ✓ 6. Infuse Stand

Made of stainless steel pipe, adjustable height with screw lock, two hanger (hook) at four different position.

#### 7. Power consumption

220 V AC / 50 Hz, ± 100 watt, back up power ± 10 minute.

#### ✓ Position

Back raise position, knee raise, Hi-Low position, Trendelenburg and reverse trendelenburg controlled with remote control.

#### ✓ Finishing

Epoxy powder coating.

#### ✓ Dimension (Approx)

Overall 2200 (L) x 990 (W) x 510 -770 (H) mm

büma